

Asan J & Image J 실무연수강좌

# Basics of Image Processing #2

박 범 우

Biomedical Imaging Infrastructure

# CONTENTS


신약개발 융합 바이오이미징 센터

Center for Bio-imaging of New Drug Development

- Image Windowing and Threshold

- Simple Transform






# CHAPTER

신약개발 융합 바이오이미징 센터  
Center for Bio-imaging of New Drug Development



## Image Windowing and Threshold

# Windowing & Threshold



Original Image



Contrast Stretched Image  
By Windowing




Binary (Mask) Image  
By Threshold

Windowing is  
**displaying technique**



Threshold is  
**selection technique**



# Contrast Stretching




# Negative Image






# Windowing




Display image **brighter** than original




Display image **darker** than original

# Brightness & Contrast





# Various Terminology for Windowing



Count: 262144

Min: 35

Mean: 128.229

Max: 240

StdDev: 42.764

Mode: 154 (3171)

- Lower than Min be displayed Black
- Upper than Max be displayed White

- **Min**


- **Max**

- **Window Level(WL)**




- **Window Width(WW)**

- **Contrast**

- **Brightness**




# Example for Windowing on ImageJ





1. Min-Max
  1. Min: 0
  2. Max: 255
2. Window Level & Width
  1. WL: 128
  2. WW: 255
3. Contrast & Brightness
  1. Contrast: 1
  2. Brightness: 0

Try to use Lena!




# Threshold



**Threshold Value** (Signal intensity level) is important role to separate object from background.



# Example for Threshold #1



cthead1.png


Original image



Threshold  
lower = 100 upper = 255

Binary image

# Example for Threshold #2




Original image





Binary image

# Example for Threshold

-Lena-



Original image




Binary image




# Example for Threshold

-blob-




Original image



Overlay Mask



Binary image








# CHAPTER

신약개발 융합 바이오이미징 센터  
Center for Bio-imaging of New Drug Development



## Simple Transform

# Transformation






# Flip

Original image




Flip  
Horizontally



Flip  
Vertically



Flip  
Vertically &  
Horizontally




# Rotate

Original image



90° Right



90° Left




Original image



15° Right  
Enlarge image



15° Right





# Translate


Original image



40 pixels  
X-axis




40 pixels  
Y-axis





# Scaling


Original image



For X-axis  
x0.5 scaling

For Y-axis  
x0.5 scaling





# Thank you